

Product Datasheet

Rat GDF5 ELISA Kit (orb439728)

Description	Rat Growth Differentiation Factor 5 (GDF5) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	M0RAY4
Sensitivity	0.066ng/mL
Entrez	252835
Expiration Date	Please enquire.