

Product Datasheet

Bovine GDH ELISA Kit (orb439730)

Description	Bovine Glutamate Dehydrogenase 1 (GDH) Sandwich ELISA Kit
Reactivity	Bovine
Range	7.812-500U/L
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P00366
Sensitivity	3.4U/L
Entrez	281785
Expiration Date	Please enquire.