

Product Datasheet

Bovine GFAP ELISA Kit (orb439738)

Description	Bovine Glial Fibrillary Acidic Protein (GFAP) Sandwich ELISA Kit
Reactivity	Bovine
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q28115
Sensitivity	10.2pg/mL
Entrez	281189
Expiration Date	Please enquire.