

Product Datasheet

Human GH2 ELISA Kit (orb439748)

Description	Human Growth Hormone 2 (GH2) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P01242
Sensitivity	29.8pg/mL
Entrez	2689
Expiration Date	Please enquire.