

Product Datasheet

Bovine Growth Hormone ELISA Kit (orb439749)

Description	Bovine Growth Hormone (GH) Sandwich ELISA Kit
Reactivity	Bovine
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	B0FPG6
Sensitivity	28pg/mL
Entrez	290904
Expiration Date	Please enquire.