

Product Datasheet

Human GIF ELISA Kit (orb439773)

Description	Human Gastric Intrinsic Factor (GIF) competitive inhibition ELISA Kit
Reactivity	Human
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P27352
Sensitivity	0.95ng/mL
Entrez	2694
Expiration Date	Please enquire.