

Product Datasheet

Human Gjb1 ELISA Kit (orb439777)

Description	Human Gap Junction Protein Beta 1 (Gjb1) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P08034
Sensitivity	0.108ng/mL
Entrez	2705
Expiration Date	Please enquire.