

Product Datasheet

Human GKN2 ELISA Kit (orb439780)

Description	Human Gastrokine 2 (GKN2) Sandwich ELISA Kit
Reactivity	Human
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q86XP6
Sensitivity	6.1pg/mL
Entrez	200504
Expiration Date	Please enquire.