

Product Datasheet

Human GLDC ELISA Kit (orb439786)

Description	Human Glycine Dehydrogenase (GLDC) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P23378
Sensitivity	0.132ng/mL
Entrez	2731
Expiration Date	Please enquire.