

Product Datasheet

Mouse GLDC ELISA Kit (orb439787)

Description	Mouse Glycine Dehydrogenase (GLDC) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q91W43
Sensitivity	0.063ng/mL
Entrez	104174
Expiration Date	Please enquire.