

Product Datasheet

Rat GLUT1 ELISA Kit (orb439810)

Description	Rat Glucose Transporter 1 (GLUT1) Sandwich ELISA Kit
Reactivity	Rat
Range	1.25-80ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P11167
Sensitivity	0.49ng/mL
Entrez	24778
Expiration Date	Please enquire.