



Product Datasheet

PROM1 Antibody (orb43983)


Description	PROM1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Prominin-1 protein (508-792AA)
Target	PROM1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O43490
Dilution Range	WB:1:500-5000, IHC-P:1:1000-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Hela whole cell...



Immunohistochemistry analysis of human e...



Immunohistochemical staining of human li...