

Product Datasheet

Human gp130 ELISA Kit (orb439833)

Description	Human Glucoprotein 130 (gp130) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P40189
Sensitivity	0.115ng/mL
Entrez	3572
Expiration Date	Please enquire.