

## Product Datasheet

### Human GP1Ba ELISA Kit (orb439836)

|                          |                                                                                             |
|--------------------------|---------------------------------------------------------------------------------------------|
| <b>Description</b>       | Human Glycoprotein Ib Alpha Polypeptide,<br>Platelet (GP1Ba) Sandwich ELISA Kit             |
| <b>Reactivity</b>        | Human                                                                                       |
| <b>Range</b>             | 1.56-100ng/mL                                                                               |
| <b>Form/Appearance</b>   | Traditional                                                                                 |
| <b>Storage</b>           | All reagents should be kept according to the<br>instructions on the packaging upon receipt. |
| <b>Note</b>              | For research use only                                                                       |
| <b>Application notes</b> | serum, plasma and other biological fluids.                                                  |
| <b>Sample Types</b>      | serum, plasma and other biological fluids.                                                  |
| <b>Uniprot ID</b>        | <b>P07359</b>                                                                               |
| <b>Sensitivity</b>       | 0.63ng/mL                                                                                   |
| <b>Entrez</b>            | <b>2811</b>                                                                                 |
| <b>Expiration Date</b>   | Please enquire.                                                                             |