

Product Datasheet

Mouse GP1Bb ELISA Kit (orb439837)

Description	Mouse Glycoprotein Ib Beta Polypeptide, Platelet (GP1Bb) Sandwich ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	P56400
Sensitivity	14.7pg/mL
Entrez	14724
Expiration Date	Please enquire.