

Product Datasheet

Human GP39 ELISA Kit (orb439839)

Description	Human Glycoprotein 39, Cartilage (GP39) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	P36222
Sensitivity	12.2pg/mL
Entrez	1116
Expiration Date	Please enquire.