

Product Datasheet

Human GP73 ELISA Kit (orb439844)

Description	Human Golgi Protein 73 (GP73) Sandwich ELISA Kit
Reactivity	Human
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q8NBJ4
Sensitivity	6.5pg/mL
Entrez	51280
Expiration Date	Please enquire.