

Product Datasheet

Human GPC2 ELISA Kit (orb439851)

Description	Human Glypican 2 (GPC2) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q8N158
Sensitivity	0.063ng/mL
Entrez	221914
Expiration Date	Please enquire.