

Product Datasheet

Human GPER ELISA Kit (orb439856)

Description	Human G Protein Coupled Estrogen Receptor 1 (GPER) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q99527
Sensitivity	0.115ng/mL
Entrez	2852
Expiration Date	Please enquire.