

Product Datasheet

Human GPLD1 ELISA Kit (orb439861)

Description	Human Glycosylphosphatidylinositol Specific Phospholipase D1 (GPLD1) Sandwich ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P80108
Sensitivity	0.34ng/mL
Entrez	2822
Expiration Date	Please enquire.