

Product Datasheet

Human GREM1 ELISA Kit (orb439883)

Description	Human Gremlin 1 (GREM1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, urine and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, urine and other biological fluids.
Uniprot ID	O60565
Sensitivity	0.064ng/mL
Entrez	26585
Expiration Date	Please enquire.