

Product Datasheet

Mouse GREM1 ELISA Kit (orb439884)

Description	Mouse Gremlin 1 (GREM1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	070326
Sensitivity	0.061ng/mL
Entrez	23892
Expiration Date	Please enquire.