

Product Datasheet

Human GRIA2 ELISA Kit (orb439888)

Description	Human Glutamate Receptor, Ionotropic, AMPA 2 (GRIA2) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P42262
Sensitivity	0.115ng/mL
Entrez	2891
Expiration Date	Please enquire.