

Product Datasheet

Human GRIA3 ELISA Kit (orb439890)

Description	Human Glutamate Receptor, Ionotropic, AMPA 3 (GRIA3) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P42263
Sensitivity	0.116ng/mL
Entrez	2892
Expiration Date	Please enquire.