

Product Datasheet

Human GRIN2A ELISA Kit (orb439895)

Description	Human Glutamate Receptor, Ionotropic, N-Methyl-D-Aspartate 2A (GRIN2A) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q12879
Sensitivity	0.103ng/mL
Entrez	2903
Expiration Date	Please enquire.