

Product Datasheet

Human GRIN2C ELISA Kit (orb439898)

Description	Human Glutamate Receptor, Ionotropic, N-Methyl-D-Aspartate 2C (GRIN2C) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q14957
Sensitivity	0.118ng/mL
Entrez	2905
Expiration Date	Please enquire.