

Product Datasheet

Rat GROg ELISA Kit (orb439910)

Description	Rat Growth Regulated Oncogene Gamma (GROg) Sandwich ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q10746
Sensitivity	33pg/mL
Entrez	171551
Expiration Date	Please enquire.