

Product Datasheet

Rat gSAP ELISA Kit (orb439914)

Description	Rat Gamma-Secretase Activating Protein (gSAP) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	D3ZPQ9
Sensitivity	0.113ng/mL
Entrez	311984
Expiration Date	Please enquire.