

Product Datasheet

Human GSTa5 ELISA Kit (orb439931)

Description	Human Glutathione S Transferase Alpha 5 (GSTa5) Sandwich ELISA Kit
Reactivity	Human
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q7RTV2
Sensitivity	0.25ng/mL
Entrez	221357
Expiration Date	Please enquire.