

Product Datasheet

Mouse GSTo1 ELISA Kit (orb439934)

Description	Mouse Glutathione S Transferase Omega 1 (GSTo1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	009131
Sensitivity	0.051ng/mL
Entrez	14873
Expiration Date	Please enquire.