

Product Datasheet

Human GSTt1 ELISA Kit (orb439937)

Description	Human Glutathione S Transferase Theta 1 (GSTt1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P30711
Sensitivity	0.055ng/mL
Entrez	2952
Expiration Date	Please enquire.