

Product Datasheet

Human GUCA2A ELISA Kit (orb439944)

Description	Human Guanylate Cyclase Activator 2A (GUCA2A) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000fmol/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q02747
Sensitivity	8.78fmol/mL
Entrez	2980
Expiration Date	Please enquire.