



Product Datasheet

CD41 antibody (orb43996)

Description	Mouse monoclonal antibody to CD41
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Leukocytes of patient suffering from LGL-type leukaemia.
Target	CD41
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 2 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P08514
Entrez	3674
Dilution Range	Flow cytometry: Recommended dilution: 2 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human CD41 positive thrombocytes.



Flow cytometry surface staining pattern ...