

Product Datasheet

Human H2AFV ELISA Kit (orb439966)

Description	Human H2A Histone Family, Member V (H2AFV) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q71UI9
Sensitivity	0.057ng/mL
Entrez	94239
Expiration Date	Please enquire.