

Product Datasheet

Human HABP2 ELISA Kit (orb439975)

Description	Human Hyaluronan Binding Protein 2 (HABP2) Sandwich ELISA Kit
Reactivity	Human
Range	78.12-5000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q14520
Sensitivity	28ng/mL
Entrez	3026
Expiration Date	Please enquire.