

Product Datasheet

Human HBg2 ELISA Kit (orb439990)

Description	Human Hemoglobin Gamma 2 (HBg2) Sandwich ELISA Kit
Reactivity	Human
Range	312.5-20000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, erythrocyte lysates and other biological fluids.
Sample Types	serum, plasma, erythrocyte lysates and other biological fluids.
Uniprot ID	P69892
Sensitivity	95.7ng/mL
Entrez	3048
Expiration Date	Please enquire.