

Product Datasheet

Mouse HCAM ELISA Kit (orb439997)

Description	Mouse Homing Associated Cell Adhesion Molecule (HCAM) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P15379
Sensitivity	29pg/mL
Entrez	12505
Expiration Date	Please enquire.