

Product Datasheet

Human HCII ELISA Kit (orb440001)

Description	Human Heparin Cofactor II (HCII) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	P05546
Sensitivity	0.067ng/mL
Entrez	3053
Expiration Date	Please enquire.