

Product Datasheet

Human HCRP1 ELISA Kit (orb440005)

Description	Human Hepatocellular Carcinoma Related Protein 1 (HCRP1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q8NEZ2
Sensitivity	0.049ng/mL
Entrez	137492
Expiration Date	Please enquire.