

Product Datasheet

Mouse HDAC6 ELISA Kit (orb440015)

Description	Mouse Histone Deacetylase 6 (HDAC6) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9Z2V5
Sensitivity	0.064ng/mL
Entrez	15185
Expiration Date	Please enquire.