

Product Datasheet

Mouse HEXb ELISA Kit (orb440031)

Description	Mouse Hexosaminidase B Beta (HEXb) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P20060
Sensitivity	0.118ng/mL
Entrez	15212
Expiration Date	Please enquire.