

Product Datasheet

Porcine HIF1a ELISA Kit (orb440038)

Description	Porcine Hypoxia Inducible Factor 1 Alpha (HIF1a) Sandwich ELISA Kit
Reactivity	Porcine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	A0T1H6
Sensitivity	5.8pg/mL
Entrez	396696
Expiration Date	Please enquire.