

Product Datasheet

Mouse HMG1 ELISA Kit (orb440052)

Description	Mouse High Mobility Group Protein 1 (HMG1) Sandwich ELISA Kit
Reactivity	Mouse
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P63158
Sensitivity	23.73pg/mL
Entrez	15289
Expiration Date	Please enquire.