

Product Datasheet

Porcine HMG1 ELISA Kit (orb440053)

Description	Porcine High Mobility Group Protein 1 (HMG1) Sandwich ELISA Kit
Reactivity	Porcine
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P12682
Sensitivity	28.3pg/mL
Entrez	445521
Expiration Date	Please enquire.