

Product Datasheet

Canine HMGA1 ELISA Kit (orb440055)

Description	Canine High Mobility Group AT Hook Protein 1 (HMGA1) Sandwich ELISA Kit
Reactivity	Canine
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q6URC2
Sensitivity	0.115ng/mL
Entrez	442946
Expiration Date	Please enquire.