

Product Datasheet

Human HMGN1 ELISA Kit (orb440060)

Description	Human High Mobility Group Nucleosome Binding Domain Protein 1 (HMGN1) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P05114
Sensitivity	0.114ng/mL
Entrez	3150
Expiration Date	Please enquire.