

Product Datasheet

Bovine Hpt ELISA Kit (orb440081)

Description	Bovine Haptoglobin (Hpt) Sandwich ELISA Kit
Reactivity	Bovine
Range	156.25-10000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q2TBU0
Sensitivity	41.5ng/mL
Entrez	280692
Expiration Date	Please enquire.