

Product Datasheet

Canine Hpt ELISA Kit (orb440082)

Description	Canine Haptoglobin (Hpt) Sandwich ELISA Kit
Reactivity	Canine
Range	62.5-4000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P19006
Sensitivity	24.5ng/mL
Expiration Date	Please enquire.