

Product Datasheet

Rat Hpt ELISA Kit (orb440086)

Description	Rat Haptoglobin (Hpt) Sandwich ELISA Kit
Reactivity	Rat
Range	15.625-1000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P06866
Sensitivity	7.3ng/mL
Entrez	24464
Expiration Date	Please enquire.