

Product Datasheet

Rat HPX ELISA Kit (orb440088)

Description	Rat Hemopexin (HPX) Sandwich ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P20059
Sensitivity	32pg/mL
Entrez	58917
Expiration Date	Please enquire.